

DISPENSING MEDICATION USING THE ADAPTA CAP & SYRINGE



- 1** Shake bottle well.
- 2** Remove outside white cap.
- 3** Insert flat surface end of plastic stopper so it is level with top of bottle.
- 4** Insert syringe tip into hole of plastic stopper.
- 5** Turn bottle over and draw prescribed amount of liquid into syringe.
- 6** Turn bottle back over and pull syringe out. Leave stopper in bottle.
NOTE: medication will not leak out of syringe
- 7** Replace white cap on bottle and dispense medication.



Stokes Pharmacy
TECHNOLOGY AND TRADITION

Please call the pharmacy if you have any questions, 800-754-5222