

Dispensing Medication Using the Adapta Cap & Syringe

- 1 Shake bottle well.
- 2 Remove outside white cap.
- 3 Insert tapered end of plastic stopper into bottle and push in until the stopper is level with top of bottle.
- 4 Insert syringe tip into hole of plastic stopper.
- 5 Turn bottle over and draw prescribed amount of liquid into syringe.
- 6 Turn bottle back over and pull syringe out. Leave stopper in bottle.
NOTE: Do not wash syringe.
- 7 Replace white cap on bottle and dispense medication.



Questions? Call 800.754.5222