

Medi-Turn™

Transdermal Dial-Up Pen

Though you should not come in direct contact with this medication, for your safety, we recommend protective gloves be worn when applying.



Priming for First Time Use

1. Remove protective cap.
2. Turn base until a small amount of medication exits the applicator tip.
3. Continue turning until the blue arrows are aligned (Don't worry about wasting medication. Extra is included for this purpose).
4. Wipe the applicator tip with a clean tissue. Re-cap the pen as you prepare for administering the dose.

Regular Dosing

1. Remove the protective cap and verify the blue arrows are properly aligned.
2. Unless otherwise directed by your veterinarian, rotate the base of the pen 2 FULL REVOLUTIONS using the aligned arrows as a marker.
3. Apply the gel dispensed directly to the inner pinna of the animal's ear.
4. Spread the medication over the entire ear surface. Be careful not to get medication into the ear canal. Alternate ears with each dose.
5. Wipe the applicator tip with a clean tissue after application and replace the cap.
6. Wash hands.
7. If you have other animals or small children, keep cat in separate area for 15 minutes while drug is absorbed.

Clean ears twice weekly for buildup. Use damp cloth or solution of alcohol and water with cotton ball to clean. Wash hands after.